



Protect Your Crew and Patients with Integrated Vehicle Intelligence.

AEV is taking ambulance safety to the next level with Traumahawk Telematics, which comes standard on many of our Type I, Type III and Medium-Duty ambulances. This intelligent system provides real-time monitoring of your ambulance fleet, giving you the vital stats you need to enhance driver and vehicle safety. From maintenance alerts to fuel use monitoring and GPS location services, Traumahawk Telematics puts you in the driver's seat when it comes to elevating the safety and efficiency of your ambulance fleet. This smart monitoring technology comes with a one-year subscription to ACETECH™ Vehicle Intelligence, giving your ambulances the ability to send real-time safety alerts to nearby motorists' navigation apps and in-vehicle infotainment systems.

Protect the lives of your ambulance crews and patients by elevating your fleet's performance with Traumahawk Telematics.

TRAUMAHAWK TELEMATICS™ FEATURES

- Real-time driver safety coaching with voice alerts
- GPS location services
- Accelerometer and turning data
- Vehicle accessory information (e.g., lights on)
- Maintenance alerts
- Battery monitoring
- Fuel use monitoring
- Ambulance status
- Geofencing
- OBD2 codes
- Collision alerts
- Data collection, preserved and stored (even in cellular dead zones)
- One-year free subscription to ACETECH™ iNSight cloud-based software
- HAAS Alert Safety Cloud® digital alerting capabilities via ACETECH CONNECT
- Powered by ACETECH



The American-Crafted Experience™

At the heart of AEV is our unwavering commitment to caring for our customers. We are inspired by a calling to build ambulances with the highest standards of American-crafted excellence and approach every relationship with equally high standards. We do that by manufacturing superior products and delivering an outstanding experience with a welcoming spirit and warmth that extends long after the sale.



AEV PRODUCT EXCELLENCE

We expertly craft the broadest range of ambulances available today, making AEV the number one selling brand in America.



QUALITY & SAFETY BUILT IN

With just-in-time, LEAN production processes and a focus on quality and safety, our talented craftspeople ensure superior value goes into every AEV ambulance.



AFTER SALES SUPPORT AND WARRANTY

At AEV, the purchase is just the beginning. We provide outstanding warranties, including up to 7-year/70,000-mile paint and electrical warranties and up to a 20-year structural warranty.

We proudly stand behind our products and care for you long after the sale with a dedicated network of independent dealers, backed by your AEV team. Our experienced team is ready to help you. Contact our sales team at **1-800-374-9749**, or find your local dealer at **www.aev.com/dealer-locator**.



Jefferson, NC 28640 | 1-800-374-9749 | www.aev.com

© 2024 REV Group, Inc. All rights reserved. | 6/2024
All trademarks are the property of their respective owners.